

Lifter

Model AWL-150

Bruno's Economy Lifts

Lifting Capacity: 200 lb/90 kg



Transport Your Scooter or Manual Wheelchair in the Trunk or Boot with a Fold-Down Lift!

Bruno's Lifter conveniently lifts and stores your Personal Mobility Device (manual folding wheelchair, travel scooter, travel powerchair) inside your vehicle with the touch of a button. Connect the docking device, push a button and let your Lifter do the heavy lifting for you!

- Works with minivans, SUVs, Crossovers, full-size vans and sedans (check with your local Bruno dealer to confirm compatibility with your vehicle).
- No modifications to your Personal Mobility Device or vehicle
- Power Raise/Lower and Manual Rotate Operation
- Can reinstall in your next vehicle
- Up to \$1200 is available to help pay for a vehicle lift through Mobility Reimbursement Programs



Lift folds down when not in use

Important Note: Confirm all applications with an authorized Bruno dealer

Represented by:

BRUNO INDEPENDENT LIVING AIDS, INC.®

1780 Executive Drive, Oconomowoc, WI 53066

www.bruno.com

