

# BetterLiving

## PressureSurface Double Layer Extension



### Dual layer 200 mm long foam extension for use with 150 mm high pressure surfaces.

Suitable for use in most standard residential care, hospital and home care beds with sleeping platform extension capabilities. High density memory foam upper layer for additional pressure relief under “at risk” heel zone.

Medium density foam core base. Double strength multi-stretch waterproof cover reduces the potential for shear injuries.

ITEM **Double Layer Extension**

CODE **BSPS0472**

CAPACITY NA

HEIGHT 150 mm

WIDTH 865 mm

LENGTH 200 mm

TOP COVER Latex free polyurethane

WARRANTY 6 year body, 2 year cover

COMPLIANCE **ARTG 176837**

CARE To clean, wipe down the top cover and PVC base with warm water and a solution of sodium hypochlorite or similar (up to 10,000 ppm available chlorine). Top cover can be machine washed up to 80°C.



200 mm extension for beds with sleeping platform extension



Multi-density foam layers for maximum strength and comfort