



COMPLEMENT YOUR QUICKIE PULSE WITH A JAY® CUSHION

Order your JAY® Cushion shipped complete on-chair with any Quickie Pulse chair.



Pulse 5™

JAY® Cushions

Comfort Cushions



JAY Soft Combi™ P



JAY Ion™



JAY Union™

High needs Cushions



JAY J2™/J2DC



JAY J3™



JAY Fusion™



JAY Easy™



Comfort and performance
in a proven design



Captain's Seat
model shown

Quickie Pulse Specifications

SPECIFICATIONS	PULSE 5 BC	PULSE 5 CC
Drive Wheel Position	Mid	Mid
Speed	8 km/ h	8 km/ h
Overall Weight	107 kg	113 kg.
Overall Width	24"	24"
Overall Length	34"	34"
User Weight Capacity	136 kg	136 kg
Seat Widths	16" - 20"	16", 18", 20"
Seat Depths	16" - 20"	15" - 22"
Front Seat To Floor Height	16" - 20"	16.25" - 18.25" (from pan)
Seating Frame Options	ASAP I	Captain's Seat
Transit Approved	ISO	WC-19
Warranty	Limited/Lifetime	Limited/Lifetime
Shroud Colours	Candy Apple Red, Candy Blue	Candy Apple Red, Candy Blue
Armrest Options	Cantilever	Cantilever
Legrest Options	Center Mount, Swing- Away	Center Mount, Swing- Away
	Manual ELR	Manual ELR
TGA #	100375	100375



© Sunrise Medical PULSE 5. 2013. Rev A.AU

Sunrise Medical Australia Pty Ltd
6 Healey Circuit, Huntingwood, NSW 2148
Phone: (02) 9678 6600
Fax: (02) 9678 6655
Email: enquiries@sunrisemedical.com.au
www.sunrisemedical.com.au

Live Without Limits

The Quickie Pulse is a compact, durable power wheelchair that offers a wide range of effective seating and electronics options.

QUICKIE® Pulse 5

Available in 2 Models:

Quickie Pulse 5 (CC): 8 Kph, 2-Pole motor, No Power Option base with Captain Seat
Quickie Pulse 5 (BC): 8 Kph, 2-Pole motor, Tilt and ASAP I, Adjustable Seat Frame



POWER BASE DESIGNED FOR PERFORMANCE

- Mid wheel drive
- Reliable and compact for indoor accessibility
- Flexible suspension system for optimal outdoor performance
- Sleek base style available with shroud colour options: Candy Apple Red or Candy Blue



ELECTRONICS OPTIONS

VR-2™ Non-Expandable Controls

- Industry standard non-expandable electronics
- Dependable, highly flexible, and simple to use
- 2 button Thru Drive for Tilt models




MULTIPLE SEAT FRAME OPTIONS Captain Seat (CC)

- Fully adjustable seat size, legrest placement and armrest width
- Convertible seat pan for flexibility in positioning solutions
- Sleek, comfortable style
- Semi-recline and fold-down backrest
- Mounting accessories
 - O₂ Holder
 - Crutch Holder
 - Walker Holder
- Available in 16", 18" and 20" width and 15"- 22" depth
- Seat to floor height adjustment from 16"- 20"



ASAP I Seat frame (BC)

ADDITIONAL BACKREST OPTIONS FOR ASAP I SEAT FRAMES (BC)



JAY J3 Backrest



SPOT Curved Backrest

ASAP SEAT (BC):

- ASAP seat frame with static tilt adjustment of 0° and 4°
- SPOT Tilt on ASAP seat frame provides up to 50° of tilt
- Backrests: Jay Curved back or Jay 3 back
- In-field width, depth and back height adjustability
16 -20" Wide, 16- 20" Deep

The JAY J3™ back's stylish design, quick-release hardware, and ability to meet a multitude of clinical needs make it an excellent complement to all Pulse models.

The SPOT Curved back's solid support system and option of Whitmyer headrests and positioning components make it a fantastic option on the Pulse.