



# XCEL

## HD MID-WHEEL DRIVE

The New Heavy Duty Mid-wheel Power Wheelchair



QUICKIE XCEL offers a standard weight capacity of 249 kg. The new Mid-wheel Heavy Duty power base features an Advance Geometric Design (A.G.D) which helps enhance the stability, mobility and comfort for the end user. As part of the Quickie X family of Canadian-made wheelchairs, the Xcel benefits through it's GC3 technology; Gravity Compensation, Ground Compensation and the Gravity Control without compromise on stability, control and freedom.



MADE IN CANADA



### Mid-wheel Drive Featuring:

- GC3 Technology
- Short Turning Radius
- PG Drives Technology
- Fully Adjustable Eclipse Seat Frame
- Low Seat to Floor Height
- Extra Duty Components





# Heavy Duty Mid-Wheel Drive Power Wheelchair

The Quickie® Xcel™ promise: Improving People's Lives. Unsurpassed quality that provides the bariatric user to better access of their environment.

## GC<sup>3</sup> TECHNOLOGY™



### 1. Gravity Control

Stability and traction without necessary adjustments to compensate for various body types relative to complex or dynamic seating.



### 2. Ground Compensation

While riding over uneven surfaces, the hard movements and vibrations are controlled and not absorbed by the seat or user.



### 3. Gravity compensation

The unique suspension will distribute the weight over the entire footprint and assure an anterior and posterior stability.

## SEATING



### Bariatric Fully Adjustable Rehab Seat

- Width adjustable 20" to 28"
- Depth adjustable 17" to 24"

#### Options:

- Back angle 84° to 120° every 6°
- HD power Tilt 45°
- HD power recline with mechanical shear reduction and power elevating & articulating legrests
- HD power ELR/ALR footboard



Complement your Quickie® Xcel™ with **JAY J3** seating

## FEATURES

Weight Capacity	249 KG Standard
Front Caster	8" Caster
Drive Wheel	14 " Cast wheel pneumatic with or without insert
Transit	4-Point wheelchair transport included
Electronics	Standard: R-net EL 90 Amp expandable controller Options: R-net PM 120 amp fully expandable controller
Seat to floor height	16.5" – with standard seating) 17 1/4 " – with tilt
Max Speed	8.5 KPH
Overall Base Width	70 CM
Overall Base Length	110.5 CM
Seating	HD fully adjustable seat bariatric Eclipse
Warranty	5 years for frame, please see owner's manual for all warranty information
Colours available	Standard: Black Options: TBA



[www.QuickieXPower.ca](http://www.QuickieXPower.ca)



© Sunrise Medical. XCEL.2013.Rev A.AU

Sunrise Medical Australia Pty Ltd  
6 Healey Circuit, Huntingwood, NSW 2148  
Phone: (02) 9678 6600  
Fax: (02) 9678 6655  
Email: [enquiries@sunrisemedical.com.au](mailto:enquiries@sunrisemedical.com.au)  
[www.sunrisemedical.com.au](http://www.sunrisemedical.com.au)