



XPERIENCE MID-WHEEL DRIVE




is ideal for an active end user sharing their time between indoor and outdoor activities. This Mid-Wheel Drive power wheelchair allows for easy maneuvering over day to day obstacles while offering unsurpassed stability provided by the GC3 suspension technology.



MADE IN CANADA 

Mid-Wheel Drive Featuring:

- GC3 Technology
- Short Turning Radius
- PG Drives Technology
- Fully Adjustable ASAP or Eclipse Seat Frame
- Low Seat to Floor Height





Mid-Wheel Drive Power Wheelchair

The Quickie® Xperience™ has stable maneuverability which revolutionizes the standards in power wheelchairs.

GC³ TECHNOLOGY™



1. Gravity Control

Stability and traction without necessary adjustments to compensate for various body types relative to complex or dynamic seating.

2. Ground Compensation

While riding over uneven surfaces, the hard movements and vibrations are controlled and not absorbed by the seat or user.

3. Gravity compensation

The unique suspension will distribute the weight over the entire footprint and assure an anterior and posterior stability.

SEATING



Eclipse Adjustable Rehab Seat

- Width adjustable 10" to 24"
- Depth adjustable 10" to 22"
- Back angle 86° to 122° every 4°
- Power Tilt 50° Tilt
- Power recline with mechanical shear reduction and power elevating & articulating legrests.
- Power ELR/ALR footboard.



Complement your Quickie® Xperience™ with JAY J3 seating

FEATURES Quickie® Xperience™

Weight Capacity	136 KG Standard, 181 KG Heavy Duty option
Front Caster	6" Caster, 8" Caster pneumatic with or without insert
Drive Wheel	14" Cast wheel pneumatic with or without insert
Transit	4-Point wheelchair transport included
Electronics	R-Net EL 90 Amp expandable R-Net 120 Amp expandable, OPTION
Seat to floor height	40.6 CM (16" *with 14" drive wheels)
Max Speed	10 KPH
Overall Base Width	64.8 CM
Overall Base Length	110.5 CM
Seating	Fully adjustable seat: Eclipse
Warranty	5 years for frame, please see owner's manual for all warranty information
Colours available	Black, Candy Red, Razzleberry, Burnt Orange, Red Pearl and Blue Pearl



www.QuickieXPower.ca



© Sunrise Medical.XPERPSE.2013.Rev A.AU

Sunrise Medical Australia Pty Ltd
6 Healey Circuit, Huntingwood NSW 2148
Phone: (02) 9678 6600
Fax: (02) 9831 2244
Email: enquiries@sunrisemedical.com.au
www.sunrisemedical.com.au