



XPLORE

REAR-WHEEL DRIVE



is an excellent choice for an active user spending time outdoors. The directional stability at high speeds in all types of environments, helps ensure that the end user's optimal Center of Gravity will be maintained without any adjustments required.



Rear-Wheel Drive Featuring:

- GC3 Technology
- Short Turning Radius
- PG Drives Technology Electronics
- Fully Adjustable ASAP or Eclipse Seat Frame
- Low Seat to Floor Height

MADE IN CANADA 





Rear-Wheel Drive Power Wheelchair

The Quickie® Xplore™ has stable maneuverability which revolutionizes the standards in power wheelchairs.

GC³ TECHNOLOGY™



1. Gravity Control

Stability and traction without necessary adjustments to compensate for various body types relative to complex or dynamic seating.



2. Ground Compensation

While riding over uneven surfaces, the hard movements and vibrations are controlled and not absorbed by the seat or user.



3. Gravity compensation

The unique suspension will distribute the weight over the entire footprint and assure an anterior and posterior stability.

SEATING



ASAP Adjustable Seat

- Width adjustable 16" to 20"
- Depth adjustable 16" to 20"

Options

- Back angle 86° to 122° every 4°
- HD Power Tilt 50°



Complement your Quickie® Xplore™ with JAY J3 seating

FEATURES Quickie® Xplore™

Weight Capacity	136 KG Standard, 181 KG Heavy Duty option
Front Caster	8" Caster, 9" Caster pneumatic with or without insert
Drive Wheel	14" Cast wheel pneumatic with or without insert
Transit	c4-Point wheelchair transport included
Electronics	R-Net EL 90 Amp fully expandable controller R-net 120 Amp fully expandable controller OPTION
Seat to floor height	16.5"
Max Speed	14 KPH
Overall Base Width	65 CM
Overall Base Length	110 CM
Seating	Fully adjustable seat: Eclipse
Warranty	5 years for frame, please see owner's manual for all warranty information
Colours available	Black, Red, Pearl, Blue Pearl and Burnt Orange



www.QuickieXPower.ca



© Sunrise Medical. XPLOPER.2013.Rev A.AU

Sunrise Medical Australia Pty Ltd
6 Healey Circuit, Huntingwood NSW 2148
Phone: (02) 9678 6600
Fax: (02) 9678 6655
Email: enquiries@sunrisemedical.com.au
www.sunrisemedical.com.au